

Anti-CD89 antibody [MIP7c] ab53357

[1 References](#) [1 图像](#)

概述

产品名称	Anti-CD89抗体[MIP7c]
描述	小鼠单克隆抗体[MIP7c] to CD89
宿主	Mouse
特异性	ab53357 recognises the EC1 domain of the CD89 extracellular region.
经测试应用	适用于: Inhibition Assay
种属反应性	与反应: Human
免疫原	Soluble CD89
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
存储溶液	Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	MIP7c
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab53357于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Inhibition Assay		

应用说明

Flow Cyt: Use 1µg/100µl as primary antibody.
 Inhib: Use 4µg for 10⁵ cells.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Binds to the Fc region of immunoglobulins alpha. Mediates several functions including cytokine production.

组织特异性

Isoform A.1, isoform A.2 and isoform A.3 are differentially expressed between blood and mucosal myeloid cells. Isoform A.1, isoform A.2 and isoform A.3 are expressed in monocytes. Isoform A.1 and isoform A.2 are expressed in alveolar macrophages; however only one isoform is expressed at alveolar macrophages surfaces.


序列相似性

Contains 2 Ig-like C2-type (immunoglobulin-like) domains.

细胞定位

Secreted and Cell membrane.

图片



Neutrophil phagocytosis of IgA immune complexes in presence or absence (control) of ab53357.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors