

Anti-CFP10 antibody ab45074

★★★★★ **3 Abreviews** **5 References**

概述

产品名称	Anti-CFP10抗体
描述	兔多克隆抗体to CFP10
宿主	Rabbit
特异性	ab45074 recognises native and E coli expressed recombinant CFP10.
经测试应用	适用于: WB
种属反应性	与反应: Mycobacterium tuberculosis
免疫原	Recombinant full length protein corresponding to CFP10.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
纯化说明	Purity: > 95% (by SDS-PAGE)
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab45074于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	Use at an assay dependent concentration.

靶标

相关性	CFP10 (10 kDa culture filtrate protein) is a low molecular weight protein isolated from Mycobacterium tuberculosis culture filtrates. CFP10 and the 6 kDa early secretory antigen of T cells (ESAT6) are secreted in abundance by Mycobacterium tuberculosis and are frequently recognized by T cells from infected people. The genes encoding these proteins have been deleted from the genome of the vaccine strain Bacille Calmette Guerin (BCG). CFP10 is one of a small number of secreted M. tuberculosis proteins now identified, that appear to play key roles in tuberculosis pathogenesis and in the stimulation of immunity.
细胞定位	Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors