



Anti-CGA/hCG alpha antibody ab196505

2 图像

概述

产品名称	Anti-CGA/hCG alpha抗体
描述	兔多克隆抗体to CGA/hCG alpha
宿主	Rabbit
经测试应用	适用于: WB, ICC/IF, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human CGA/hCG alpha aa 1 to the C-terminus. Database link: P01215
	 Run BLAST with  Run BLAST with
阳性对照	JAR cell lysate, HepG2 cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 49% PBS, 0.87% Sodium chloride, 50% Glycerol (glycerin, glycerine)</p> <p>PBS is without Mg^{2+} and Ca^{2+}</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab196505于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/2000. Predicted molecular weight: 13 kDa.
ICC/IF		1/50 - 1/200.
IHC-P		Use at an assay dependent concentration.

靶标


序列相似性

Belongs to the glycoprotein hormones subunit alpha family.

细胞定位


Secreted.

图片



Anti-CGA/hCG alpha antibody (ab196505) at 1/500 dilution + JAR cell lysate

Predicted band size: 13 kDa



Immunofluorescent analysis of HepG2 cells labeling hCG alpha with ab196505 at 1/50 dilution. Lower image shows DAPI for nuclear staining (blue).

Immunocytochemistry/ Immunofluorescence - Anti-CGA/hCG alpha antibody (ab196505)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors