

# Anti-CGA/hCG alpha antibody [INN-hFSH-132] ab11232

★★★★★ [1 Abreviews](#) [9 References](#)

## 概述

产品名称	Anti-CGA/hCG alpha抗体[INN-hFSH-132]
描述	小鼠单克隆抗体[INN-hFSH-132] to CGA/hCG alpha
宿主	Mouse
经测试应用	适用于: ELISA, IP, WB, RIA, IHC-P
种属反应性	与反应: Human
免疫原	Chemical/ Small Molecule corresponding to Human CGA/hCG alpha. Human Follicle Stimulating Hormone Database link: <a href="#">P01215</a>
表位	ab11232 recognises the alpha 4 epitope of the alpha-subunit of Human Chorionic Gonadotrophin (hCG), hLH, hFSH, hTSH and free alpha-subunit.
常规说明	<p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	INN-hFSH-132

骨髓瘤	x63-Ag8.653
同种型	IgG1

## 应用

**The Abpromise guarantee** Abpromise™ 承诺保证使用ab11232于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		
IP		
WB		
RIA		
IHC-P	★★★★★ (1)	

## 应用说明

ELISA: Use at an assay dependent dilution.  
 IHC-P: Used at a dilution of 1/1500 for 2 hrs on human placenta sections (see Abreview for further details).  
 IP: Use at an assay dependent dilution.  
 RIA: Use at an assay dependent dilution.  
 WB: Use at an assay dependent dilution (see Berger et al).  
  
 Not tested in other applications.  
 Optimal dilutions/concentrations should be determined by the end user.

## 靶标

序列相似性	Belongs to the glycoprotein hormones subunit alpha family.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors