

### Anti-Chk1 (phospho S286) antibody ab62243

★★★★★ [2 Abreviews](#) [1 图像](#)

#### 概述

产品名称	Anti-Chk1 (phospho S286)抗体
描述	兔多克隆抗体to Chk1 (phospho S286)
宿主	Rabbit
经测试应用	适用于: ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human Chk1 aa 250-350 (phospho S286). Database link: <a href="#">O14757</a>
阳性对照	HeLa cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride</p>
纯度	Without Mg <sup>2+</sup> and Ca <sup>2+</sup>
纯化说明	Immunogen affinity purified
克隆	ab62243 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site.
同种型	多克隆
	IgG

## 应用


**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab62243于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/1000.

## 靶标

## 图片



ab62243 at 1/500 - 1/1000 dilution staining Chk1 in HeLa cells by Immunofluorescence, in the absence or presence of the immunising peptide.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors