

# Anti-Chromogranin A antibody [LK2H10] - BSA and Azide free ab212290

## 2 图像

### 概述

产品名称	Anti-Chromogranin A抗体[LK2H10] - BSA and Azide free
描述	小鼠单克隆抗体[LK2H10] to Chromogranin A - BSA and Azide free
宿主	Mouse
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 不与反应: Sheep, Rabbit, Guinea pig
免疫原	Tissue, cells or virus corresponding to Human Chromogranin A. (Human pheochromocytoma).
阳性对照	Human adrenal gland and small cell lung carcinoma tissues; PC12 cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Constituent: 100% PBS
无载体	是
纯度	Protein A/G purified
纯化说明	ab212290 was purified from Bioreactor Concentrate by Protein A/G.
克隆	单克隆
克隆编号	LK2H10

同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab212290于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 0.25 - 0.5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. (Primary incubation for 30 minutes at RT).

靶标


功能	Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.
序列相似性	Belongs to the chromogranin/secretogranin protein family.
翻译后修饰	Sulfated on tyrosine residues and/or contains sulfated glycans. O-glycosylated with core 1 or possibly core 8 glycans.
细胞定位	Secreted. Neuroendocrine and endocrine secretory granules.
形式	According to the Swiss-Prot database, this protein has many different forms which correspond to various molecular weights; thus the possible variation from the predicted molecular weight of 51kDa.

图片



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human adrenal gland tissue labeling Chromogranin A with ab212290 at 0.5 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Chromogranin A antibody [LK2H10] - BSA and Azide free (ab212290)



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human small cell lung carcinoma tissue labeling Chromogranin A with ab212290 at 0.5 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Chromogranin A antibody [LK2H10] - BSA and Azide free (ab212290)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors