

Anti-CNTF antibody ab28951

[1 References](#) [1 图像](#)

概述

产品名称	Anti-CNTF抗体
描述	鸡多克隆抗体to CNTF
宿主	Chicken
经测试应用	适用于: WB
种属反应性	与反应: Mouse
免疫原	Two peptides for human CNTF were mixed and used as immunogens. Both peptides are highly conserved against human and mouse sequences.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.02% Sodium azide Constituent: PBS
纯化说明	This chicken polyclonal antibody is a crude IgY preparation from egg yolks using a precipitation method (PEG, Polyethylene Glycol).
克隆	多克隆
同种型	IgY

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab28951于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/2000. Predicted molecular weight: 24 kDa.

靶标

功能	CNTF is a survival factor for various neuronal cell types. Seems to prevent the degeneration of motor axons after axotomy.
组织特异性	Nervous system.
序列相似性	Belongs to the CNTF family.
细胞定位	Cytoplasm.

图片



Anti-CNTF antibody (ab28951) at 1/1000 dilution + Mouse brain extract at 40ug total protein.

Predicted band size: 24 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors