

Anti-COX7B antibody [2G7H8AD9] ab110269

1 References **1 图像**

概述

产品名称	Anti-COX7B抗体[2G7H8AD9]
描述	小鼠单克隆抗体[2G7H8AD9] to COX7B
宿主	Mouse
特异性	ab110269 cross-reacts weakly with Human and Rat.
经测试应用	适用于: WB
种属反应性	与反应: Cow
免疫原	Full length protein. This information is considered to be commercially sensitive.
阳性对照	Bovine heart lysate
常规说明	<p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> <p>Product was previously marketed under the MitoSciences sub-brand.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.5</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: HEPES buffered saline</p>
纯化说明	Near homogeneity as judged by SDS-PAGE. ab110269 was produced in vitro using hybridomas grown in serum-free medium, and then purified by biochemical fractionation
克隆	单克隆

克隆编号	2G7H8AD9
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab110269于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 25 µg/ml. Predicted molecular weight: 9 kDa.

靶标

功能	This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.
序列相似性	Belongs to the cytochrome c oxidase VIIb family.
细胞定位	Mitochondrion inner membrane.

图片



Anti-COX7B antibody [2G7H8AD9] (ab110269) at 25 µg/ml + Bovine heart lysate at 1 µg

Predicted band size: 9 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors