

Anti-CPNE7 antibody [CPNE7-01] ab68547

1 图像

概述	
产品名称	Anti-CPNE7抗体[CPNE7-01]
描述	小鼠单克隆抗体[CPNE7-01] to CPNE7
宿主	Mouse
特异性	ab68547 recognizes C terminus of human CPNE7.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment (GST-tag) corresponding to Human CPNE7 (C terminal). Database link: Q9UBL6
阳性对照	WB: HeLa cell lysate- nuclear fraction
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.098% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	CPNE7-01
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™** 承诺保证使用ab68547于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 70 kDa (predicted molecular weight: 70 kDa).

靶标

功能	May function in membrane trafficking. Exhibits calcium-dependent phospholipid binding properties.
序列相似性	Belongs to the copine family. Contains 2 C2 domains. Contains 1 VWFA domain.

图片



All lanes : Anti-CPNE7 antibody [CPNE7-01] (ab68547) at 1 µg/ml

Lane 1 : HeLa cell lysate- nuclear fraction

Lane 2 : HeLa cell lysate- cytoplasmic fraction

Predicted band size: 70 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors