

Anti-CRMP5 antibody ab36203

★★★★★ [1 Abreviews](#) [6 References](#) [2 图像](#)

概述

产品名称	Anti-CRMP5抗体
描述	兔多克隆抗体to CRMP5
宿主	Rabbit
特异性	No cross reactivity was detected with CRMP1, CRMP2, CRMP3 or CRMP4.
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	A 15 residue synthetic peptide derived from C-terminal domain of human collapsing-response mediator.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab36203于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/2000. Predicted molecular weight: 65 kDa.
IHC-P		Use at an assay dependent concentration.

靶标

图片



All lanes : Anti-CRMP5 antibody (ab36203)

Lane 1 : E.coli lysate transformed with cDNA of human CRMP1

Lane 2 : E.coli lysate transformed with cDNA of human CRMP2


Lane 3 : E.coli lysate transformed with cDNA of human CRMP3

Lane 4 : E.coli lysate transformed with cDNA of human CRMP4

Lane 5 : E.coli lysate transformed with cDNA of human CRMP5

Predicted band size: 65 kDa

Observed band size: 65 kDa



Ab36203 staining human CRMP5 in olfactory receptor neurons within the olfactory mucosa by immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors