

Anti-CYFIP2 antibody ab95969

★★★★☆ [2 Abreviews](#) [9 References](#) [2 图像](#)

概述

产品名称	Anti-CYFIP2抗体
描述	兔多克隆抗体to CYFIP2
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to a region within amino acids 346-610 of Human CYFIP2 (UniProt ID: Q96F07 isoform 2; NP_055191).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 79.99% PBS, 20% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab95969于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB	★★★★★ (1)	1/500 - 1/3000. Predicted molecular weight: 146 kDa.
IHC-P		1/100 - 1/500.

靶标

相关性 CYFIP2 is involved in T cell adhesion and p53 dependent induction of apoptosis. It does not bind RNA but is up regulated significantly in CD4+ T lymphocytes from patients with multiple sclerosis. There are 2 isoforms produced by alternative splicing.

细胞定位 Cytoplasm perinuclear region. Cell junction, synapse, synaptosome. Treatment with leptomycin-B triggers translocation to the nucleus.


图片



Anti-CYFIP2 antibody (ab95969) at 1/1000 dilution + MOLT4 whole cell lysate at 30 µg

Predicted band size: 146 kDa

5% SDS PAGE



ab95969, at a 1/100 dilution, staining CYFIP2 in paraffin embedded NCIN87 xenograft by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CYFIP2 antibody (ab95969)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors