

# Anti-DAP1 antibody [E59] ab32056

**重组** RabMAb

## 1 References [4 图像](#)

### 概述

产品名称	Anti-DAP1抗体[E59]
描述	兔单克隆抗体[E59] to DAP1
宿主	Rabbit
经测试应用	<b>适用于:</b> WB, IHC-P, ICC/IF <b>不适用于:</b> Flow Cyt
种属反应性	<b>与反应:</b> Mouse, Human
免疫原	Synthetic peptide within Human DAP1 aa 50-150 (C terminal). The exact sequence is proprietary.
阳性对照	NIH3T3 cell lysate, Hela cell lysate and human stomach carcinoma
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"><li>- High batch-to-batch consistency and reproducibility</li><li>- Improved sensitivity and specificity</li><li>- Long-term security of supply</li><li>- Animal-free production</li></ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> <p>Rat: We have preliminary internal testing data to indicate this antibody may not react with this species. Please contact us for more information.</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 49% PBS, 50% Glycerol (glycerin, glycerine), 0.05% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	E59

同种型

IgG

## 应用

**The Abpromise guarantee** Abpromise™ 承诺保证使用 ab32056 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/10000. Detects a band of approximately 11 kDa.
IHC-P		1/100 - 1/250. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		1/250.

## 应用说明


Is unsuitable for Flow Cyt.

## 靶标

## 功能


Involved in mediating interferon-gamma-induced cell death.

## 图片



Immunofluorescent analysis of DAP1 expression in HeLa cells using 1/250 ab32056.

Immunocytochemistry/ Immunofluorescence - Anti-DAP1 antibody [E59] (ab32056)




Anti-DAP1 antibody [E59] (ab32056) at 1/10000 dilution + NIH3T3 cell lysate

**Observed band size:** 11 kDa

**Additional bands at:** 23 kDa. We are unsure as to the identity of these extra bands.

Western blot - Anti-DAP1 antibody [E59] (ab32056)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DAP1 antibody [E59] (ab32056)

ab32056, at a 1/100 dilution, staining human stomach carcinoma by immunohistochemistry, Paraffin embedded tissue.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

### Why choose a recombinant antibody?



Anti-DAP1 antibody [E59] (ab32056)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors