


### Anti-DAP12 antibody ab93846

★★★★★ [1 Abreviews](#) [1 References](#) [2 图像](#)

#### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-DAP12抗体  |
| 描述    | 兔多克隆抗体to DAP12  |
| 宿主    | Rabbit  |
| 经测试应用 | 适用于: WB<br>不适用于: ICC/IF   |
| 种属反应性 | 与反应: Human<br>预测可用于: Sheep, Cow, Dog   |
| 免疫原   | Synthetic peptide corresponding to Human DAP12 aa 50 to the C-terminus conjugated to keyhole limpet haemocyanin.<br>(Peptide available as <a href="#">ab107906</a> )  |
| 阳性对照  | This antibody gave a positive signal in the following lysates: Lung (Human) Tissue Lysate.  |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.   |
| 存储溶液 | pH: 7.40<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS   |
|      | Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help. |

## 应用

**Abpromise™** 承诺保证使用ab93846于以下的经测试应用


| 应用 | Ab评论 | 说明   |
|----|------|--|
| WB |      | Use a concentration of 1 µg/ml. Detects a band of approximately 12 kDa (predicted molecular weight: 12 kDa). |

应用说明 Is unsuitable for ICC/IF.

靶标

|       |   |
|-------|---|
| 功能    | Non-covalently associates with activating receptors of the CD300 family. Cross-linking of CD300-TYROBP complexes results in cellular activation.  |
| 组织特异性 | Expressed at low levels in the early development of the hematopoietic system and in the promonocytic stage and at high levels in mature monocytes. Expressed in hematological cells and tissues such as peripheral blood leukocytes and spleen. Also found in bone marrow, lymph nodes, placenta, lung and liver. Expressed at lower levels in different parts of the brain especially in the basal ganglia and corpus callosum.  |
| 疾病相关  | Defects in TYROBP are a cause of polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (PLOSL) [MIM:221770]; also called presenile dementia with bone cysts or Nasu-Hakola disease (NHD). PLOSL is a recessively inherited disease characterized by a combination of psychotic symptoms rapidly progressing to presenile dementia and bone cysts restricted to wrists and ankles. PLOSL has a global distribution, although most of the patients have been diagnosed in Finland and Japan, with an estimated population prevalence of 2x10 <sup>-6</sup> in the Finns. |
| 序列相似性 | Belongs to the TYROBP family.   |
| 翻译后修饰 | Tyrosine phosphorylated.  |
| 细胞定位  | Membrane.   |

图片



**All lanes :** Anti-DAP12 antibody (ab93846) at 1 µg/ml

**Lane 1 :** Wild-type THP-1 cell lysate

**Lane 2 :** TYROBP CRISPR-Cas9 edited THP-1 cell lysate


Lysates/proteins at 40 µg per lane.

Performed under reducing conditions.

**Predicted band size:** 12 kDa

**Observed band size:** 12 kDa

False colour image of Western blot: Anti-DAP12 antibody staining at 1 µg/ml, shown in black; Mouse anti-GAPDH antibody [6C5] ([ab8245](#)) loading control staining at 1/20000 dilution, shown in red. In Western blot, ab93846 was shown to bind specifically to DAP12. A band was observed at 12 kDa in wild-type THP-1 cell lysates with no signal observed at this size in TYROBP CRISPR-Cas9 edited cell line [ab273726](#) (CRISPR-Cas9 edited cell lysate [ab273774](#)). The band observed in the CRISPR-Cas9 edited lysate lane below 12 kDa is likely to represent a truncated form of DAP12. This has not been investigated further and the functional properties of the gene product have not been determined. To generate this image, wild-type and TYROBP CRISPR-Cas9 edited THP-1 cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 5 % BSA in TBS-0.1 % Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4 °C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times before development with Optiblot (ECL reagent [ab133456](#)) and imaged with 2 minutes exposure time. Secondary antibodies used were HRP conjugated Goat anti-Rabbit (H+L) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed ([ab216776](#)) at 1/20000 dilution.



Western blot - Anti-DAP12 antibody (ab93846)

Anti-DAP12 antibody (ab93846) at 1 µg/ml + Lung (Human) Tissue Lysate at 10 µg

### Secondary

Goat polyclonal Secondary Antibody to Rabbit IgG - H&L (HRP), pre-adsorbed at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 12 kDa

**Observed band size:** 12 kDa

**Additional bands at:** 200 kDa. We are unsure as to the identity of these extra bands.

**Exposure time:** 12 minutes

The DAP12 protein sequence contains a 21-residue signal sequence at the N-terminus. Cleavage of this peptide may account for the doublet observed.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors