

Anti-DCBLD2/ESDN antibody ab224102

2 图像

概述

产品名称	Anti-DCBLD2/ESDN抗体
描述	兔多克隆抗体to DCBLD2/ESDN
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Recombinant fragment corresponding to Human DCBLD2/ESDN aa 550-700. Database link: Q96PD2 <div style="text-align: right;">  Run BLAST with  Run BLAST with </div>
阳性对照	IHC-P: Human heart tissue. ICC/IF: U-2 OS cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab224102 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml. Fixation/Permeabilization: PFA/Triton X-100.

靶标

组织特异性

Highly expressed in testis, heart, skeletal muscle and also in cultured vascular smooth muscle cells.


序列相似性

Contains 1 CUB domain.
Contains 1 F5/8 type C domain.
Contains 1 LCCL domain.

细胞定位


Membrane.

图片



PFA-fixed, Triton X-100 permeabilized U-2 OS (human bone osteosarcoma epithelial cell line) cells stained for DCBLD2/ESDN (green) using ab224102 (4 µg/ml) in ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-DCBLD2/ESDN antibody (ab224102)



Paraffin-embedded human heart tissue stained for DCBLD2/ESDN with ab224102 (1/20) in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DCBLD2/ESDN antibody (ab224102)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors