

Anti-DEFB-4 antibody [L13-10-D1] ab14419

★★★★★ [1 Abreviews](#) [4 References](#)

概述

产品名称	Anti-DEFB-4抗体[L13-10-D1]
描述	小鼠单克隆抗体[L13-10-D1] to DEFB-4
宿主	Mouse
经测试应用	适用于: ELISA, IHC-P, RIA, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human DEFB-4 aa 25-61. Sequence: ELDRICGYGTARCRKKCRSQEYRIGRCPNTYACCLRK) Database link: Q8WTQ1 (Peptide available as ab69499)  Run BLAST with  Run BLAST with
阳性对照	IHC-P: Uterus and kidney sections.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 0.79% Tris HCl
纯度	Protein G purified
克隆	单克隆
克隆编号	L13-10-D1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab14419于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration.
RIA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Can be blocked with <u>Recombinant human DEFB-4 protein (ab69499)</u> .

靶标

功能	Has antimicrobial activity. Synergistic effects with lysozyme and DEFB103.
组织特异性	High expression in the testis. Gastric antrum exhibited relatively high levels. A lower expression is observed in uterus and neutrophils thyroid gland, lung, and kidney. No detectable expression in other tissues tested.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors