

Anti-Dengue Virus 1+2+3+4 antibody ab26837

★★★★★ [1 Abreviews](#) [3 References](#) [1 图像](#)

概述

产品名称	抗登革病毒1+2+3+4抗体
描述	兔多克隆抗体to登革病毒1+2+3+4
宿主	Rabbit
经测试应用	适用于: ICC/IF
种属反应性	与反应: Dengue virus 1, Dengue virus 2, Dengue virus 3, Dengue virus 4
免疫原	Viral particles from C6/36 [mosquito cell line] amplified dengue type 2 virus.
阳性对照	dengue type III infected BHK21 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
纯度	Whole antiserum
纯化说明	Ammonium sulphate precipitated
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab26837于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF	★★★★★ (1)	

应用说明

ICC/IF: Use at a dilution of 1/50-1/100.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Dengue and dengue hemorrhagic fever (DHF) are caused by one of four closely related, but antigenically distinct, virus serotypes (DEN1, DEN2, DEN3, and DEN4), of the genus Flavivirus.

图片



Immunofluorescent analysis of dengue type 3 infected BHK21 cells, using 10 µg/ml ab26837. Note that not every cells are infected.

Immunocytochemistry/ Immunofluorescence - Anti-Dengue Virus 1+2+3+4 antibody (ab26837)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors