

# Anti-Dengue Virus NS1 glycoprotein antibody [DN1] ab41490

[1 References](#) [2 图像](#)

## 概述

产品名称	Anti-Dengue Virus NS1 glycoprotein抗体[DN1]
描述	小鼠单克隆抗体[DN1] to Dengue Virus NS1 glycoprotein
宿主	Mouse
特异性	Recognises the non-structural protein NS1 from Dengue Virus 2 (16681 strain). Reacts only weakly with NS1 from Dengue Virus 2 NGC strain, use <a href="#">ab41623</a> for detection of Dengue Virus 2 NGC NS1.
经测试应用	适用于: WB
种属反应性	与反应: Dengue virus 2
免疫原	Full length native protein purified from Dengue Virus 2 infected supernatant
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Proclin 150 Constituents: 10% BSA, 89.9% RPMI 1640
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	DN1
同种型	IgG2a
轻链类型	kappa

应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab41490于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		1/50 - 1/100. Use under non reducing condition. Detects a band of approximately 40 kDa (predicted molecular weight: 40 kDa).


靶标

**相关性**

NS1 is one of 7 Dengue Virus non-structural proteins which are thought to be involved in viral replication. NS1 exists as a monomer in its immature form but is rapidly processed in the endoplasmic reticulum to form a stable dimer. A small amount of NS1 remains associated with intracellular organelles where it is thought to be involved in viral replication. The rest of NS1 is found either associated with the plasma membrane or secreted as a soluble hexadimer. NS1 is essential for viral viability but its precise biological function is unknown. Antibodies raised in response to NS1 in viral infection can cross react with cell surface antigens on epithelial cells and platelets and this has been implicated in the development of Dengue Hemorrhagic fever.

图片





**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors