



Anti-DFFB/CAD antibody ab8408

2 图像

概述	
产品名称	Anti-DFFB/CAD抗体
描述	兔多克隆抗体to DFFB/CAD
宿主	Rabbit
经测试应用	适用于: ICC/IF, WB
种属反应性	与反应: Human
免疫原	<p>Synthetic peptide corresponding to Human DFFB/CAD aa 3-18.</p> <p>Sequence:</p> <p>QKPKSVKLRALRSPKRC</p> <p>Database link: O76075 (Peptide available as ab8440)</p> <div> Run BLAST with  Run BLAST with</div>
阳性对照	K562 or Jurkat whole cell lysate A 40 kDa band can be detected.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p>
纯度	Affinity purified
克隆	多克隆
同种型	IgG
轻链类型	unknown

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab8408于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use a concentration of 5 µg/ml.
WB		1/500 - 1/2000. Detects a band of approximately 40 kDa.

靶标

功能	Nuclease that induces DNA fragmentation and chromatin condensation during apoptosis. Degrades naked DNA and induces apoptotic morphology.
序列相似性	Contains 1 CIDE-N domain.
细胞定位	Cytoplasm. Nucleus.

图片




All lanes : Anti-DFFB/CAD antibody (ab8408) at 1/500 dilution

Lane 1 : Jurkat whole cell lysate

Lane 2 : K562 whole cell lysate

Observed band size: 40 kDa



ab8408 at 5µg/ml staining DFFB/CAD in Jurkat cells by ICC/IF.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors