

Anti-Dlx5 antibody ab64827

★★★★☆ [4 Abreviews](#) [4 References](#) [2 图像](#)

概述

产品名称	Anti-Dlx5抗体
描述	兔多克隆抗体to Dlx5
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from an internal sequence of human Dlx5.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 50% Glycerol, 0.87% Sodium chloride</p>
纯度	Without Mg ²⁺ and Ca ²⁺
纯化说明	Immunogen affinity purified
克隆	ab64827 was affinity purified from rabbit antiserum by affinity chromatography using epitope specific immunogen.
同种型	多克隆
	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab64827于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.
WB		1/500 - 1/1000. Detects a band of approximately 32 kDa (predicted molecular weight: 32 kDa).

靶标

细胞定位

Nuclear

图片



All lanes : Anti-Dlx5 antibody (ab64827) at 1/500 dilution

Lane 1 : HeLa cell extract at 30 µg


Lane 2 : 293 cell extract at 30 µg

Lane 3 : COLO205 cell extract at 30 µg

Lane 4 : COLO205 cell extract at 30 µg
plus immunising peptide at 10 µg

Predicted band size: 32 kDa

Observed band size: 32 kDa



ab64827, at 1/50 dilution, staining Dlx5 in paraffin embedded human lung carcinoma tissue by Immunohistochemistry in the absence or presence of the immunising peptide.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors