

Anti-DNaseI antibody [4B10] ab138710

★★★★★ [1 Abreviews](#) [1 图像](#)

概述

产品名称	Anti-DNaseI抗体[4B10]
描述	小鼠单克隆抗体[4B10] to DNaseI
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Cow, Human
免疫原	Recombinant full length Cow DNase I; produced using CHO-based proprietary suspension cell line; purified from cell culture supernatant.
阳性对照	Recombinant Cow DNase I expressed into the supernatant of CHO cell culture medium.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 100% PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	4B10
同种型	IgG1

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用ab138710于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000. Predicted molecular weight: 31 kDa.

图片



Anti-DNaseI antibody [4B10] (ab138710) at 1/5000 dilution +
Recombinant Cow DNase I expressed into the supernatant of CHO
cell culture medium at 10 µl

Predicted band size: 31 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors