

# Anti-DNTTIP1/TDIF1 antibody ab174663

[1 References](#) [2 图像](#)


### 概述

产品名称	Anti-DNTTIP1/TDIF1抗体
描述	兔多克隆抗体to DNTTIP1/TDIF1
宿主	Rabbit
经测试应用	适用于: IP, WB
种属反应性	与反应: Human 预测可用于: Cow, Orangutan 
免疫原	Synthetic peptide within Human DNTTIP1/TDIF1 aa 279-329. The exact sequence is proprietary. NP_443183.1 Sequence: ASDDYRGCLDLKLEELKSFVLPSWMVEKMRKYMETLRTENEH RAVEAPPQ T  Database link: <a href="#">Q9H147</a> <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
阳性对照	Whole cell lysate from HeLa, 293T and Jurkat cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate





Western blot - Anti-DNTTIP1/TDIF1 antibody  
(ab174663)

**All lanes** : Anti-DNTTIP1/TDIF1 antibody (ab174663) at  
0.100000001490116 µg/ml

**Lane 1** : HeLa whole cell lysate

**Lane 2** : 293T whole cell lysate

**Lane 3** : Jurkat whole cell lysate

Lysates/proteins at 50 µg per lane.

Developed using the ECL technique.

**Predicted band size:** 37 kDa

**Exposure time:** 3 minutes

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors