



DyLight® 550 Anti-DDDDK tag (Binds to FLAG® tag sequence) antibody [M2] ab117495

概述

产品名称	DyLight® 550荧光Anti-DDDDK tag (Binds to FLAG® tag sequence)抗体[M2]
描述	DyLight® 550荧光小鼠单克隆抗体[M2] to DDDDK tag (Binds to FLAG® tag sequence)
宿主	Mouse
偶联物	DyLight® 550. Ex: 562nm, Em: 576nm
经测试应用	适用于: ICC/IF
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to DDDDK tag (Binds to FLAG® tag sequence). <div>  Run BLAST with ExPASy  Run BLAST with NCBI </div>
常规说明	<p>FLAG® is a registered trade mark of Sigma Aldrich Biotechnology LP. It is used here for informational purposes only.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 1.6% Sodium phosphate, 0.88% Sodium chloride, 1% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	M2
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab117495于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/500 - 1/5000.

靶标

相关性

This is a useful tool for the localisation and characterisation of DDDDK tagged proteins (Binds to FLAG® tag sequence).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors