

### Anti-E tag antibody ab95868

★★★★★ [1 Abreviews](#) [2 References](#) [1 图像](#)

#### 概述

产品名称	Anti-E tag抗体
描述	山羊多克隆抗体to E tag
宿主	Goat
经测试应用	适用于: WB
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to E tag conjugated to keyhole limpet haemocyanin. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	WB: E. coli whole cell lysate expressing a multi-tag fusion protein.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 6.8 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab95868于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent concentration. For colorimetric techniques use at a dilution of 1/1000 - 1/10000 and for chemiluminescent techniques use at a dilution of 1/1000 - 1/30000.

## 靶标

## 相关性

E tag is a peptide epitope tag, sequence GAPVPYPDPLEPR.

## 图片



**All lanes :** Anti-E tag antibody (ab95868) at 0.04 µg/ml

**Lane 1 :** E. coli whole cell lysate expressing a multi-tag fusion protein at 0.2 µg

**Lane 2 :** E. coli whole cell lysate expressing a multi-tag fusion protein at 0.1 µg

**Lane 3 :** E. coli whole cell lysate expressing a multi-tag fusion protein at 0.05 µg

**Exposure time:** 30 seconds

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors