

Anti-EBNA 3B antibody ab16127

1 References

概述

产品名称	Anti-EBNA 3B抗体
描述	羊多克隆抗体to EBNA 3B
宿主	Sheep
特异性	This antibody is specific for the Epstein Barr Virus Nuclear Antigen 3B.
经测试应用	适用于: ICC, ELISA, IHC-P, Dot blot, WB
种属反应性	与反应: Epstein-Barr virus
免疫原	Recombinant full length protein corresponding to Epstein-Barr virus EBNA 3B.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.08% Sodium azide Constituent: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	Purified by ammonium sulfate precipitation This antibody is provided as a 0.2µm sterile filtered solution.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab16127于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 165 kDa.

靶标

相关性 EBNA 3B is one of only 10 viral proteins expressed in Epstein Barr Virus (EBV) transformed B lymphoblasts cultured in vitro, though the expression of EBNA 3B is not required for EBV-mediated transformation in vitro. This 165kDa nuclear protein functions as a transcription factor.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors