

# Anti-EBV gp340/220 Envelope Protein antibody [10B5] - BSA and Azide free ab6525

## 2 References

### 概述

产品名称	Anti-EBV gp340/220 Envelope蛋白抗体[10B5] - BSA and Azide free
描述	小鼠单克隆抗体[10B5] to EBV gp340/220 Envelope蛋白- BSA and Azide free
宿主	Mouse
特异性	This antibody is specific for Epstein-Barr Virus.
经测试应用	适用于: ELISA, ICC/IF
种属反应性	与反应: Epstein-Barr virus
免疫原	Tissue, cells or virus corresponding to EBV gp340/220 Envelope Protein.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 100% PBS
无载体	是
纯度	Protein A purified
克隆	单克隆
克隆编号	10B5
骨髓瘤	NS1/1-Ag4-1
同种型	IgG1

## 应用

**The Abpromise guarantee** Abpromise™承诺保证使用ab6525于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/800.
ICC/IF		Use at an assay dependent dilution.

## 靶标

**相关性** Epstein-Barr virus is a member of the herpesvirus family and one of the most common human viruses. The envelope glycoprotein Gp340/Gp220 is the most abundant vomponent of the viral envelope and is believed to be responsible for EBV binding to CR2 receptor on human B-Cells.

**细胞定位** Virion membrane: most abundant component of the viral envelope.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors