

# Anti-ELL3 antibody ab67415

**3 References**   [1 图像](#)

### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-ELL3抗体   |
| 描述    | 小鼠多克隆抗体to ELL3  |
| 宿主    | Mouse   |
| 特异性   | Reactive against mammalian transfected lysate.  |
| 经测试应用 | 适用于: WB   |
| 种属反应性 | 与反应: Human  |
| 免疫原   | Recombinant full length protein within Human ELL3 aa 1-400. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please <b><a href="#">contact</a></b> our Scientific Support team to discuss your requirements.  |
| 阳性对照  | Transfected 293T lysate   |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | pH: 7.4<br>Constituent: PBS   |
| 纯度   | Protein A purified  |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

## 应用

### The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab67415 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度; 实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab 评论 | 说明   |
|----|-------|--|
| WB |       | 1/500 - 1/1000. Detects a band of approximately 60 kDa (predicted molecular weight: 45 kDa). |

## 靶标


### 相关性

ELL3 is a testis-specific RNA polymerase II elongation factor which increases the catalytic rate of transcription elongation. ELL3 is approximately 50% similar to both ELL and ELL2.

### 细胞定位

Nuclear

## 图片



Western blot - Anti-ELL3 antibody (ab67415)

**All lanes :** Anti-ELL3 antibody (ab67415) at 1/500 dilution

**Lane 1 :** ELL3 transfected 293T lysate

**Lane 2 :** Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size:** 45 kDa

**Observed band size:** 43.78 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors