

Anti-ESAT6 antibody ab45073

★☆☆☆☆ **3 Abreviews** **4 References**

概述

产品名称	Anti-ESAT6抗体
描述	兔多克隆抗体to ESAT6
宿主	Rabbit
特异性	The rabbit polyclonal antibody to EsaT-6 recognizes the EsaT-6 (Rv3875) protein of Mycobacterium tuberculosis and Mycobacterium bovis.
经测试应用	适用于: WB
种属反应性	与反应: Mycobacterium bovis, Mycobacterium tuberculosis
免疫原	Recombinant full length protein corresponding to ESAT6.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab45073于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★☆☆☆☆ (2)	Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 11 kDa. Use under reducing conditions.

靶标

相关性	The ESAT6 antigen from Mycobacterium tuberculosis is a dominant target for cell mediated immunity in the early phase of tuberculosis (TB) in TB patients as well as in various animal models. ESAT6 is not found in M. bovis BCG. Its function is not known but it elicits high level of IFN-gamma from memory effector cells during the first phase of a protective immune response.
细胞定位	Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors