

Anti-FAK (phospho Y397) antibody [EP2160Y] ab81298

重组 RabMAb

★★★★★ **2 Abreviews** **118 References** **8 图像**

概述

产品名称	Anti-FAK (phospho Y397)抗体[EP2160Y]
描述	兔单克隆抗体[EP2160Y] to FAK (phospho Y397)
宿主	Rabbit
特异性	The activation of phosphorylation of FAK is reported to be related to developmental processes in brain tissue (PMID: 14642275, PMID: 21118706). So expression level of phosphorylated modified FAK in normal brain is quite low causing not easy to be detected. We suggest optimizing experimental protocols (increasing lysate amount, using lower dilution or higher sensitivity ECL substrate) to improve results.
经测试应用	适用于: WB, ICC/IF 不适用于: IHC-P
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: Human, mouse and rat brain tissue, treated NIH/3T3 with 10mM Pervanadate for 60 min and treated HeLa with 10mM Pervanadate for 60 min cell lysates. ICC/IF: SK-N-SH cell line.
常规说明	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA

纯度	Protein A purified
克隆	单克隆
克隆编号	EP2160Y
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab81298于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB	★★★★★ (1)	1/1000. Predicted molecular weight: 119 kDa.
ICC/IF	★★★★★ (1)	Use a concentration of 5 µg/ml.

应用说明 Is unsuitable for IHC-P.

靶标

功能	Non-receptor protein-tyrosine kinase implicated in signaling pathways involved in cell motility, proliferation and apoptosis. Activated by tyrosine-phosphorylation in response to either integrin clustering induced by cell adhesion or antibody cross-linking, or via G-protein coupled receptor (GPCR) occupancy by ligands such as bombesin or lysophosphatidic acid, or via LDL receptor occupancy. Microtubule-induced dephosphorylation at Tyr-397 is crucial for the induction of focal adhesion disassembly. Plays a potential role in oncogenic transformations resulting in increased kinase activity.
组织特异性	Expressed in all organs tested, in lymphoid cell lines, but most abundantly in brain.
序列相似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. FAK subfamily. Contains 1 FERM domain. Contains 1 protein kinase domain.
结构域	The first Pro-rich domain interacts with the SH3 domain of CRK-associated substrate (BCAR1) and CASL. The carboxy-terminal region is the site of focal adhesion targeting (FAT) sequence which mediates the localization of FAK1 to focal adhesions.
翻译后修饰	Phosphorylated on 6 tyrosine residues upon activation. Microtubule-induced dephosphorylation at Tyr-397 could be catalyzed by PTPN11 and regulated by ZFYVE21. Dephosphorylated by PTPN11 upon EPHA2 activation by its ligand EFNA1.
细胞定位	Cell junction > focal adhesion. Cell membrane. Constituent of focal adhesions.

图片



All lanes : Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298) at 1/1000 dilution

Lane 1 : Untreated NIH/3T3 (Mouse embryonic fibroblast) whole cell lysate

Lane 2 : NIH/3T3 treated with 10mM Pervanadate for 60 min whole cell lysate

Lane 3 : Untreated HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysate

Lane 4 : HeLa treated with 10mM Pervanadate for 60 min whole cell lysate

Lane 5 : Mouse brain tissue lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 119 kDa

Observed band size: 119 kDa


Exposure time: 40 seconds

Blocking buffer and concentration: 5% NFDm/TBST.

Diluting buffer and concentration: 5% NFDm/TBST.

The activation of phosphorylation of FAK is reported to be related to developmental processes in brain tissue (PMID: 14642275, PMID: 21118706). So expression level of phosphorylated modified FAK in normal brain is quite low causing not easy to be detected. We suggest optimizing experimental protocols (increasing lysate amount, using lower dilution or higher sensitivity ECL substrate) to improve results.

[ab181602](#) has been used as a loading control.



All lanes : Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298) at 1/2000 dilution (purified)

Lane 1 : Untreated mouse brain

Lane 2 : Mouse brain treated with alkaline phosphatase

Lysates/proteins at 10 µg per lane.

Secondary


All lanes : HRP conjugated goat anti-rabbit IgG (H+L) at 1/1000 dilution

Predicted band size: 119 kDa


Observed band size: 119 kDa

Blocking Buffer: 5% NFDM/TBST

Dilution Buffer: 5% NFDM/TBST



ICC/IF image of unpurified ab81298 stained SK-N-SH (human neuroblastoma cell line) cells. The cells were 100% methanol fixed (5 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab81298, 10µg/ml) overnight at +4°C in PBS containing 1% BSA and 0.1% tween. The secondary antibody (green) was **ab96899**, DyLight® 488 goat anti-rabbit IgG (H+L) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.



All lanes : Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298) at 1/1000 dilution (purified)

Lane 1 : Untreated rat brain

Lane 2 : Rat brain treated with alkaline phosphatase

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP conjugated goat anti-rabbit IgG (H+L) at 1000 µg

Predicted band size: 119 kDa

Observed band size: 119 kDa


Blocking Buffer: 5% NFDM/TBST

Dilution Buffer: 5% NFDM/TBST

Immunocytochemistry/ Immunofluorescence - Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298)

Unpurified ab81298 staining FAK (phospho Y397) in SK-N-SH (human neuroblastoma cell line) cells treated with anandamide (in water soluble emulsion) ([ab120429](#)), by ICC/IF. Increase in FAK (phospho Y397) expression correlates with increased concentration of anandamide (in water soluble emulsion), as described in literature.

The cells were incubated at 37°C for 5 minutes in media containing different concentrations of [ab120429](#) (anandamide (in water soluble emulsion)) in water, fixed with 100% methanol for 5 minutes at -20°C and blocked with PBS containing 10% goat serum, 0.3 M glycine, 1% BSA and 0.1% tween for 2h at room temperature. Staining of the treated cells with ab81298 (5 µg/ml) was performed overnight at 4°C in PBS containing 1% BSA and 0.1% tween. A DyLight 488 goat anti-rabbit polyclonal antibody ([ab96899](#)) at 1/250 dilution was used as the secondary antibody. Nuclei were counterstained with DAPI and are shown in blue. Membranes are stained in red with WGA.



All lanes : Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298) at 1/2000 dilution (unpurified)

Lane 1 : human brain tissue lysates, untreated.

Lane 2 : human brain tissue lysates treated with AP.

Lane 3 : rat brain tissue lysates, untreated.

Lane 4 : rat brain tissue lysates treated with AP.

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP labelled goat anti rabbit. at 1/2000 dilution

Predicted band size: 119 kDa

Observed band size: 119 kDa

Immunocytochemistry/ Immunofluorescence - Anti-FAK (phospho Y397) antibody [EP2160Y] (ab81298)

Unpurified ab81298 staining FAK (phospho Y397) in SK-N-SH (human neuroblastoma cell line) cells treated with anandamide (ethanol solution) ([ab120087](#)), by ICC/IF. Increase in FAK (phospho Y397) expression correlates with increased concentration of anandamide (ethanol solution), as described in literature. The cells were incubated at 37°C for 10 minutes in media containing different concentrations of [ab120087](#) (anandamide (ethanol solution)) in ethanol, fixed with 100% methanol for 5 minutes at -20°C and blocked with PBS containing 10% goat serum, 0.3 M glycine, 1% BSA and 0.1% tween for 2h at room temperature. Staining of the treated cells with ab81298 (5 µg/ml) was performed overnight at 4°C in PBS containing 1% BSA and 0.1% tween. A DyLight 488 goat anti-rabbit polyclonal antibody ([ab96899](#)) at 1/250 dilution was used as the secondary antibody. Nuclei were counterstained with DAPI and are shown in blue. Membranes are stained in red with WGA.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-FAK (phospho Y397) antibody [EP2160Y]
(ab81298)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors