

Anti-FAM129B antibody ab75638

1 图像

概述

产品名称	Anti-FAM129B抗体
描述	兔多克隆抗体to FAM129B
宿主	Rabbit
经测试应用	适用于: ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic non-phosphopeptide derived from human FAM129B around the phosphorylation site of serine 712 (QVSS ^P P).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, PBS</p>
纯度	Without Mg ²⁺ and Ca ²⁺
纯化说明	Immunogen affinity purified
克隆	ab75638 was purified using an epitope-specific immunogen
同种型	多克隆
	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab75638于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF		1/100 - 1/200.

靶标

序列相似性

Belongs to the Niban family.
Contains 1 PH domain.

图片



ab75638 at 1/100 dilution staining FAM129B in human HeLa cells by ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-FAM129B antibody (ab75638)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors