

Anti-FAM13A antibody ab122440

[3 References](#) [5 图像](#)

概述

产品名称	Anti-FAM13A抗体
描述	兔多克隆抗体to FAM13A
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to Human FAM13A aa 904-973 (C terminal). Sequence: TDFSARCFLDQFEDDADGFISPMDDKIPSKCSQDTGLSNLHA ASIPELLE HLQEMREEKKRIRKLRDFE Database link: O94988  Run BLAST with  Run BLAST with
阳性对照	IHC-P: Human testis, small intestine, Fallopian tube and cerebral cortex Tissue ICC/IF: Human cell line U-2 OS
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab122440于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/500 - 1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.

靶标


组织特异性

Isoform 1 is widely expressed, with highest expression in skeletal muscle, thymus, brain and lung. Isoform 3 is less abundant than isoform 1 and predominantly expressed in kidney, pancreas, liver, lung and thymus.

序列相似性


Belongs to the FAM13 family.
Contains 1 Rho-GAP domain.

图片




Immunohistochemical analysis of human small intestine labeling FAM13A with ab122440 showing strong cytoplasmic positivity in glandular cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM13A antibody (ab122440)




Immunocytochemistry/ Immunofluorescence - Anti-FAM13A antibody (ab122440)

Immunocytochemistry/ Immunofluorescence analysis of U2-OS cells labeling FAM13A with ab122440 showing localization to nucleoli, cytosol & cell junctions.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM13A antibody (ab122440)

Immunohistochemical analysis of human fallopian tube labeling FAM13A with ab122440 showing strong cytoplasmic positivity in glandular cells.



Immunohistochemical analysis of human cerebral cortex labeling FAM13A with ab122440 showing strong cytoplasmic positivity in neurons.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM13A antibody (ab122440)



Immunohistochemical analysis of human testis labeling FAM13A with ab122440 showing strong cytoplasmic positivity in glandular cells in seminiferous duct.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM13A antibody (ab122440)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors