

Anti-FAT4 antibody ab130076

[5 References](#) [1 图像](#)

概述

产品名称	Anti-FAT4抗体
描述	兔多克隆抗体to FAT4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Recombinant fragment
免疫原	Synthetic peptide conjugated to KLH corresponding to a region within Human FAT4 (SwissProt Q6V0I7).
阳性对照	recombinant FAT4 protein
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7 Preservative: 0.05% Sodium azide Constituents: 74% Tris buffered saline, 25% Glycerol (glycerin, glycerine), 0.5% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab130076 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 48 kDa (predicted molecular weight: 543 kDa).

靶标

功能	May function in the regulation of planar cell polarity (By similarity). Cadherins are cell-cell interaction molecules.
组织特异性	Widely expressed. Expressed in fetal brain, infant brain, brain tumor and colorectal cancer.
序列相似性	Contains 34 cadherin domains. Contains 6 EGF-like domains. Contains 2 laminin G-like domains.
细胞定位	Membrane. In the kidney, localizes to primary cilia.

图片



Anti-FAT4 antibody (ab130076) at 1/500 dilution + recombinant FAT4 protein

Predicted band size: 543 kDa

actual MW: 48kDa

Western blot - Anti-FAT4 antibody (ab130076)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors