

Anti-FCN1/M-Ficolin antibody ab70814

1 图像

概述

产品名称	Anti-FCN1/M-Ficolin抗体
描述	兔多克隆抗体to FCN1/M-Ficolin
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	3T3 cell lysate and rat kidney tissue lysate
常规说明	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3852 Ficolin-1 Antibody. 3852-100 is the same size as the 100 µg size of ab70814.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: 0.5% BSA, 30% Glycerol (glycerin, glycerine)</p>
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab70814 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab 评论	说明
WB		Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 33 , 37 kDa (predicted molecular weight: 35 kDa).

靶标

功能	Involved in serum exerting lectin activity. Binds GlcNAc.
组织特异性	Peripheral blood leukocytes. Also detected in spleen, lung, and thymus, may be due to the presence of tissue macrophages or trapped blood in these tissues. Not detected on lymphocytes and granulocytes.
序列相似性	Belongs to the ficolin lectin family. Contains 1 collagen-like domain. Contains 1 fibrinogen C-terminal domain.
细胞定位	Secreted. Found on the monocyte surface.

图片



All lanes : Anti-FCN1/M-Ficolin antibody (ab70814) at 0.5 µg/ml

Lane 1 : 3T3 cell lysate

Lane 2 : Rat kidney tissue lysate

Predicted band size: 35 kDa

Observed band size: 33,37 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors