



Anti-FGF4 antibody ab65974

6 References **2 图像**

概述

产品名称	Anti-FGF4抗体
描述	兔多克隆抗体to FGF4
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human FGF4 aa 84-97 (internal sequence). Sequence: RRL YCNVGIGFHLQ
	Database link: P08620
	 Run BLAST with  Run BLAST with
阳性对照	Hela cell lysis and human lung cancer tissue sections.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:0.2ml of distilled water will yield a concentration of 500µg/ml.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservatives: 0.025% Thimerosal (merthiolate), 0.025% Sodium azide Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab65974于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 22 kDa (predicted molecular weight: 22 kDa).
IHC-P		Use a concentration of 1 - 2 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

靶标


功能	Can transform NIH 3T3 cells from a human stomach tumor (hst) and from karposi's sarcoma (KS3). It has a mitogenic activity.
序列相似性	Belongs to the heparin-binding growth factors family.
细胞定位	Secreted.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FGF4 antibody (ab65974)

ab65974 at 1µg/ml staining FGF4 in human lung cancer tissue sections.



Anti-FGF4 antibody (ab65974) at 1 µg/ml + Hela cell lysis at 30 µg

Predicted band size: 22 kDa

Observed band size: 22 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors