

Anti-FGL2/Prothrombinase antibody ab198029

★ ★ ★ ★ ★ [1 Abreviews](#) [4 References](#) [2 图像](#)

概述

产品名称	Anti-FGL2/Prothrombinase抗体
描述	兔多克隆抗体to FGL2/Prothrombinase
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Mouse, Human
免疫原	Recombinant fragment (GST-tag) corresponding to Human FGL2/Prothrombinase (internal sequence). BC033820 Database link: Q14314
阳性对照	Mouse thymus tissue lysate;Human liver cancer tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.4 Preservative: 0.05% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 49% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab198029于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/20 - 1/100.
WB		1/450 - 1/2000. Predicted molecular weight: 50 kDa.

靶标

功能	May play a role in physiologic lymphocyte functions at mucosal sites.
组织特异性	Constitutively expressed in cytotoxic T-cells.
序列相似性	Contains 1 fibrinogen C-terminal domain.
细胞定位	Secreted.

图片




Anti-FGL2/Prothrombinase antibody (ab198029) at 1/450 dilution +
Mouse thymus tissue lysate at 40 µg

Predicted band size: 50 kDa

Exposure time: 40 seconds

8% SDS-PAGE



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue labeling FGL2/Prothrombinase with ab198029 at 1/20 dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors