

FITC Anti-6X His tag® antibody ab1206

★★★★★ [7 Abreviews](#) [28 References](#) [2 图像](#)

概述

产品名称	FITC荧光Anti-6X His tag®抗体
描述	FITC荧光兔多克隆抗体to 6X His tag®
宿主	Rabbit
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: ELISA, Flow Cyt
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to 6X His tag® conjugated to keyhole limpet haemocyanin. Run BLAST with ExPASy Run BLAST with NCBI
常规说明	Molar F/P ratio is 3.7:1. HIS-TAG® is a trademark of EMD Biosciences, Inc. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	Antibodies were immunoaffinity purified using the peptide conjugated to a solid-phase support and conjugated to fluorescein isothiocyanate.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab1206于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ELISA	★★★★★ (1)	Use at an assay dependent concentration.
Flow Cyt	★★★★★ (5)	Use at an assay dependent concentration. ab37406 - Rabbit polyclonal IgG, is suitable for use as an isotype control with this antibody.

靶标

相关性 The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. His-tags can be fused to the amino- or carboxy- termini of proteins in transfected or transformed cells.


细胞定位 Depends upon the localization of the parent protein tagged with hexahistidine.

图片



Muc1-Bi-1 and Muc1-Bi-2 can recognize Muc1 positive tumor cells.

a). Flow cytometry analysis of Muc1-Bi-1 and Muc1-Bi-2 on Muc1-negative cell lines HePG2 and CHO and Muc1-positive cell lines, HT29. Red line indicates cells with no staining; green line indicates cells with only anti-His-FITC staining (ab1206); blue line indicates cells with positive control anti-Muc1 antibody; purple line indicates cells with Muc1-Bi-1 protein; orange line indicates cells with Muc1-Bi-2 protein.



ab1206 staining 6X His tag[®] in human endometrial carcinoma cells by Flow Cytometry.

Cells were prepared by short trypsinisation in PBS/1% FBS/0.2% sodium azide and fixed with formaldehyde. The gate included all single cells. The sample was incubated with a 6X His tag[®] primary antibody for 1 hour at 4°C. ab1206 was used as the secondary. On the image, the red line represents cells stained with ab1206 alone, and the blue line is the tested sample carrying a 6XHis tag[®]. The small peak to the right represents cells that bound the tested sample.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors