

FITC Anti-6X His tag® antibody ab3554

★★★★★ [1 Abreviews](#) [1 References](#)

概述

产品名称	FITC荧光Anti-6X His tag®抗体
描述	FITC荧光鸡多克隆抗体to 6X His tag®
宿主	Chicken
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Species independent
免疫原	Synthetic peptide corresponding to 6X His tag® conjugated to keyhole limpet haemocyanin. Run BLAST with Expasy Run BLAST with NCBI
常规说明	Molar F/P ratio is 3.6:1. HIS-TAG® is a trademark of EMD Biosciences, Inc. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.00 Preservative: 0.1% Sodium azide Constituent: 0.2% BSA
纯度	Immunogen affinity purified
克隆	多克隆
骨髓瘤	unknown
同种型	IgY

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab3554于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Prior to conjugation, the antibodies were analyzed by Western Blot.
ELISA		Use at an assay dependent concentration. Prior to conjugation, the antibodies were analyzed by ELISA.

靶标

相关性

The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. His-tags can be fused to the amino- or carboxy- termini of proteins in transfected or transformed cells.

细胞定位

Depends upon the localization of the parent protein tagged with hexahistidine.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors