

FITC Anti-CD134 / OX40L receptor antibody [MRC OX-40] ab33960

1 图像

概述

产品名称	FITC荧光Anti-CD134 / OX40L receptor抗体[MRC OX-40]
描述	FITC荧光小鼠单克隆抗体[MRC OX-40] to CD134 / OX40L receptor
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	The antigen is a glycoprotein of approximately 50kD molecular weight which is related in sequence to the low affinity Nerve Growth Factor Receptor.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Rat
免疫原	Tissue, cells or virus corresponding to Rat CD134/ OX40L receptor. Phytohemagglutinin (PHA) activated Rat Lymph Node Cells.
阳性对照	Rat spleen cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	MRC OX-40

骨髓瘤 NS0/1 Ag4.1
同种型 IgG2b

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab33960于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 10µl for 10 ⁶ cells. Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl. ab91368 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.


靶标

功能 Receptor for TNFSF4/OX40L/GP34.

序列相似性 Contains 4 TNFR-Cys repeats.

细胞定位 Membrane.

图片



Staining of stimulated rat spleen cells with ab33960 by Flow Cytometry/FACS.

Flow Cytometry - FITC Anti-CD134 / OX40L receptor antibody [MRC OX-40] (ab33960)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors