


FITC Anti-CD28 antibody [CD28.2] ab134472

1 图像

概述

产品名称	FITC荧光Anti-CD28抗体[CD28.2]
描述	FITC荧光小鼠单克隆抗体[CD28.2] to CD28
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human 预测可用于: Chimpanzee, Baboon, Cynomolgus monkey, Rhesus monkey 
免疫原	DC28.1.3.3 hybridoma cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze. Store undiluted.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.14% Monobasic dihydrogen sodium phosphate, 0.87% Sodium chloride, 0.1% Gelatin
纯度	Affinity purified
克隆	单克隆
克隆编号	CD28.2
同种型	IgG1

The Abpromise guarantee

Abpromise™承诺保证使用ab134472于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. ab91356 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标

功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

图片



Flow Cytometry analysis of human peripheral blood lymphocytes labeling CD28 with ab134472 at uL (0.5 µg) (solid line) or Mouse IgG1 (FITC) isotype control at 0.5 µg (dashed line).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors