

### FITC Anti-CD28 antibody [PV-1] ab24844

#### 1 图像

#### 概述

产品名称	FITC荧光Anti-CD28抗体[PV-1]
描述	FITC荧光亚美尼亚仓鼠单克隆抗体[PV-1] to CD28
宿主	Armenian hamster
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse, Human
免疫原	Tissue, cells or virus corresponding to Mouse CD28. C57BL/6N mouse T cell lymphoma EL4 cell line
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Affinity purified
克隆	单克隆
克隆编号	PV-1
同种型	IgG

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab24844 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度; 实际最佳的稀释度/浓度应由使用者检定。

应用	Ab 评论	说明
Flow Cyt		Use 3µg for 10 <sup>6</sup> cells. <b>ab18473</b> - Armenian Hamster monoclonal IgG, is suitable for use as an isotype control with this antibody.

## 靶标

功能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
组织特异性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Membrane.

## 图片



ab24844 at a concentration of 3 µg/10<sup>6</sup> cells staining human A 20 (B cell lymphoma cells) and EL-4 cells by flow cytometry. A 50:50 mixture of the A20 B-cell lymphoma and the EL-4 T-cell lymphoma was stained with ab24844. Large cells were gated and analyzed by flow cytometry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors