

FITC Anti-CD90 / Thy1 antibody [F15-42-1] ab11155

15 References **1 图像**

概述

产品名称	FITC荧光Anti-CD90 / Thy1抗体[F15-42-1]
描述	FITC荧光小鼠单克隆抗体[F15-42-1] to CD90 / Thy1
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human CD90/ Thy1.
阳性对照	MOLT4 human T-cell lines. human T lymphoma cell line HUT78.
常规说明	<p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 1% BSA</p>
纯度	Protein G purified
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	F15-42-1

骨髓瘤

NS1

同种型

IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab11155于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Use neat. Use 10µl of the suggested working dilution to label 106 cells in 100µl.

靶标

功能

May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.


序列相似性

Contains 1 Ig-like V-type (immunoglobulin-like) domain.

细胞定位

Cell membrane.

图片



ab11155 staining CD90/Thy1 in human HUT78 T cells by Flow Cytometry analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors