

FITC Anti-Integrin alpha V beta 3 antibody [23C6] ab93513

[4 References](#) [1 图像](#)

概述

产品名称	FITC荧光Anti-Integrin alpha V beta 3抗体[23C6]
描述	FITC荧光小鼠单克隆抗体[23C6] to Integrin alpha V beta 3
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Integrin alpha V beta 3. Cell suspension containing osteoclasts from osteoclastomas.
阳性对照	M21 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS, 0.1% Gelatin
纯度	Protein G purified
克隆	单克隆
克隆编号	23C6
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab93513于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

应用说明

Flow Cyt: Use 5µl/test in 100µl volume.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


相关性

Function: The alpha-V (ITGAV) integrins are receptors for vitronectin, cytotactin, fibronectin, fibrinogen, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin and vWF. They recognize the sequence R-G-D in a wide array of ligands. ITGAV:ITGB3 binds to fractalkine (CX3CL1) and may act as its coreceptor in CX3CR1-dependent fractalkine signaling (PubMed:23125415). ITGAV:ITGB3 binds to NRG1 (via EGF domain) and this binding is essential for NRG1-ERBB signaling (PubMed:20682778). ITGAV:ITGB3 binds to FGF1 and this binding is essential for FGF1 signaling (PubMed:18441324). ITGAV:ITGB3 binds to IGF1 and this binding is essential for IGF1 signaling (PubMed:19578119). ITGAV:ITGB3 binds to PLA2G2A via a site (site 2) which is distinct from the classical ligand-binding site (site 1) and this induces integrin conformational changes and enhanced ligand binding to site 1 (PubMed:18635536, PubMed:25398877). ITGAV:ITGB3 and ITGAV:ITGB6 act as a receptor for fibrillin-1 (FBN1) and mediate R-G-D-dependent cell adhesion to FBN1 (PubMed:12807887, PubMed:17158881). (Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C (PubMed:20615244). Integrin ITGAV:ITGB5 and ITGAV:ITGB3 act as receptors for coxsackievirus A9 and B1 (PubMed:9426447, PubMed:15194773, PubMed:7519807). Integrin ITGAV:ITGB3 acts as a receptor for herpes virus 8/HHV-8 (PubMed:18045938). Integrin ITGAV:ITGB6 acts as a receptor for herpes simplex 1/HHV-1 (PubMed:24367260). Integrin ITGAV:ITGB3 acts as a receptor for Human parechovirus 1 (PubMed:11160695). Integrin ITGAV:ITGB3 acts as a receptor for West Nile virus (PubMed:23658209). In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions (PubMed:10397733). Similarity: Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats.

细胞定位

Cell Membrane. Cytoplasmic. Integrin complex. Focal adhesion.

图片



Flow Cytometry - FITC Anti-Integrin alpha V beta 3 antibody [23C6] (ab93513)

ab93513 at 5µl/test staining Integrin alpha V beta 3 in M21 cells (blue histogram) or FITC mouse IgG1, *k* Isotype control (green histogram). Total viable cells were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors