

FITC Anti-Klebsiella spp antibody ab69467

概述

产品名称	FITC荧光Anti-Klebsiella spp抗体
描述	FITC荧光兔多克隆抗体to Klebsiella spp
宿主	Rabbit
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	ab69467 reacts with Klebsiella sp. in bacterial and infected tissue samples.
经测试应用	适用于: IHC-Fr, ICC/IF
种属反应性	与反应: Klebsiella pneumoniae
免疫原	Tissue, cells or virus corresponding to Klebsiella spp. Whole cells Klebsiella pneumoniae
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium azide Constituents: 1% BSA, PBS
纯度	Affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab69467于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		1/10 - 1/50. Fix with Acetone.
ICC/IF		1/10.

靶标

相关性 Klebsiella is a genus of non-motile, Gram-negative, Oxidase-negative bacteria with a prominent polysaccharide-based capsule. Frequent human pathogens, Klebsiella organisms can lead to a wide range of disease states, notably pneumonia, urinary tract infections, septicemia, Ankylosing spondylitis, and soft tissue infections. Klebsiella pneumoniae is clinically the most important member of the Klebsiella genus of Enterobacteriaceae. New antibiotic resistant strains of K. pneumoniae are appearing, and it is increasingly found as a nosocomial infection.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors