

FITC Mouse IgG2a, monoclonal - Isotype control ab81197

4 References

概述

产品名称	FITC荧光小鼠IgG2a,单克隆抗体-同型对照
偶联物	FITC. Ex: 493nm, Em: 528nm
共轭说明	Excitation: 488 nm. Peak emission: 525 nm.
经测试应用	适用于: Flow Cyt
常规说明	<p>Reagents should be brought to room temperature ($22 \pm 3^{\circ}\text{C}$) before use. Protect cells from light sources during incubation with antibodies and prior to analysis.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C . Store at $+4^{\circ}\text{C}$ short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
存储溶液	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: 0.4% BSA, 1.74% Sodium chloride, 0.328% Sodium phosphate
纯度	IgG fraction
克隆	单克隆
克隆编号	TIB13
同种型	IgG2a

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab81197于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors