

FITC Mouse IgG2b, kappa monoclonal [MG2b-57] - Isotype control ab18427

1 References

概述

产品名称 FITC荧光小鼠IgG2b, kappa单克隆抗体[MG2b-57] -同型对照

偶联物 FITC. Ex: 493nm, Em: 528nm

经测试应用 **适用于:** Flow Cyt

免疫原 Trinitrophenl + KLH

常规说明

ab18427 MG2b-57 immunoglobulin is useful as an isotype-matched control in Western blotting, immunoprecipitation, immunohistochemical staining and immunofluorescent staining (surface or intracellular) for flow cytometric analysis.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式 Liquid

存放说明 Shipped at 4°C. Store at +4°C.

存储溶液 pH: 7.20
Preservative: 0.09% Sodium azide
Constituents: PBS, 0.2% BSA

纯度 Immunogen affinity purified

纯化说明 ab18427 antibody was purified by affinity chromatography, and conjugated with FITC under optimal conditions. The solution is free of unconjugated FITC.

Isotype control说明 ab18427 MG2b-57 immunoglobulin is useful as an isotype-matched control in Western blotting, immunoprecipitation, immunohistochemical staining and immunofluorescent staining (surface or intracellular) for flow cytometric analysis.

克隆	单克隆
克隆编号	MG2b-57
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab18427于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 20µl for 10 ⁶ cells. Can also be used in intracellular flow cytometric analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors