

FITC Mouse IgG2b, kappa [MPC-11] - isotype control ab18419

1 References

概述

产品名称 FITC荧光小鼠IgG2b, kappa [MPC-11] -同型对照

偶联物 FITC. Ex: 493nm, Em: 528nm

经测试应用 **适用于:** Flow Cyt

常规说明

This clone has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式 Liquid

存放说明 Shipped at 4°C. Store at +4°C. Store In the Dark.

存储溶液 pH: 7.20
Preservative: 0.09% Sodium azide
Constituents: PBS, 0.2% BSA

Isotype control说明 This clone has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues.

克隆 单克隆

克隆编号 MPC-11

同种型 IgG2b

轻链类型 kappa

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab18419于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Suitable for surface or intracellular staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors