

FITC Mouse IgG2b [PLPV219] - Isotype Control ab91368

1 References

概述	
产品名称	FITC荧光小鼠IgG2b [PLPV219] -同型对照
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	ab91368 was produced against a synthetic hapten, which is normally not present in Humans or animals. Therefore, ab91368 does not react with any antigen of Human or animal cells.
经测试应用	适用于: Flow Cyt
免疫原	Vegetal protein
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	PLPV219
同种型	IgG2b
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab91368于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

应用说明

Flow Cyt: Use 20µl for 10⁶ cells.
ab91368 is a negative isotype control

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors