


### FITC Anti-NCAM1 antibody [MEM-188], prediluted ab69763

#### 1 图像

#### 概述

产品名称	FITC荧光Anti-NCAM1抗体[MEM-188], prediluted
描述	FITC荧光小鼠单克隆抗体[MEM-188] to NCAM1, prediluted
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	Reacts with a 180 kDa isoform of NCAM expressed in leukocytes. Reactivity with other isoforms has not been tested.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human 预测可用于: Non human primates 
免疫原	Tissue, cells or virus corresponding to Human NCAM1. KG-1 Human acute myelogenous leukemia cell line
阳性对照	Human peripheral blood leukocytest.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.4 Preservative: 0.097% Sodium azide Constituents: 0.2% BSA, PBS
纯度	Size exclusion
纯化说明	The purified antibody is conjugated with Fluorescein isothiocyanate (FITC) under optimum conditions from >95% pure antibody (SDS PAGE). The reagent is free of unconjugated FITC and

adjusted for direct use.

克隆	单克隆
克隆编号	MEM-188
同种型	IgG2a

## 应用

**The Abpromise guarantee** **Abpromise™** 承诺保证使用ab69763于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

**应用说明** Flow Cyt: Use 20µl for 10<sup>6</sup> cells or 100µl whole blood.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.


## 靶标

**功能** This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.

**序列相似性** Contains 2 fibronectin type-III domains.  
Contains 5 Ig-like C2-type (immunoglobulin-like) domains.

**细胞定位** Secreted and Cell membrane.

## 图片



Flow Cytometry - FITC Anti-NCAM1 antibody [MEM-188], prediluted (ab69763)

ab69763 at 20µl/10<sup>6</sup> cells staining NCAM in human peripheral blood cells by flow cytometry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors