

# Anti-Flt3 / CD135 antibody [RM0021-12B6] ab73019

[2 References](#) [1 图像](#)

### 概述

产品名称	Anti-Flt3 / CD135抗体[RM0021-12B6]
描述	大鼠单克隆抗体[RM0021-12B6] to Flt3 / CD135
宿主	Rat
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse
免疫原	produced from a hybridoma (mouse myeloma) fused with spleen cells from a rat immunized with purified recombinant mouse Flt3/ CD135 extracellular domain.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Lyophilized:Reconstitute the antibody with 500 µl sterile PBS.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant. Filter sterilised.
克隆	单克隆
克隆编号	RM0021-12B6
同种型	IgG2

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab73019于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/50 - 1/100.

## 靶标

<b>功能</b>	Receptor for the FL cytokine. Has a tyrosine-protein kinase activity.
<b>组织特异性</b>	Bone marrow cells.
<b>序列相似性</b>	Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1/PDGF receptor subfamily. Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 protein kinase domain.
<b>细胞定位</b>	Membrane.

## 图片



ab73019 at 1/100 dilution, staining Flt3 / CD135 in LPS treated mouse tissue section by Immunohistochemistry (Formalin/PFA fixed paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Flt3 / CD135 antibody [RM0021-12B6] (ab73019)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors