

### Anti-FOXN3 antibody [CHESC9H4] ab50756

★★★★☆ [2 Abreviews](#) [3 References](#)

#### 概述

产品名称	Anti-FOXN3抗体[CHESC9H4]
描述	小鼠单克隆抗体[CHESC9H4] to FOXN3
宿主	Mouse
经测试应用	适用于: Dot blot, WB
种属反应性	与反应: Recombinant fragment 预测可用于: Human 
免疫原	Recombinant fragment corresponding to FOXN3.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1% BSA, PBS
纯度	Protein G purified
纯化说明	Purified using protein G column chromatography from culture supernatant of hybridoma cultured in a meium containing bovine IgG-depleted (approximately 95%) fetal bovine serum. Filtered through a 0.22 micrometer membrane.
克隆	单克隆
克隆编号	CHESC9H4
同种型	IgG1

## 应用

### The Abpromise guarantee

**Abpromise™**承诺保证使用ab50756于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Dot blot		Use at an assay dependent concentration.
WB	★★★★★ (2)	Use at an assay dependent concentration.

## 靶标

功能	Acts as a transcriptional repressor. May be involved in DNA damage-inducible cell cycle arrests (checkpoints).
序列相似性	Contains 1 fork-head DNA-binding domain.
细胞定位	Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors