

Anti-FSH beta antibody [INN-hFSH-60] ab11069

2 References

概述

产品名称	Anti-FSH beta抗体[INN-hFSH-60]
描述	小鼠单克隆抗体[INN-hFSH-60] to FSH beta
宿主	Mouse
特异性	Ab11069 Shows strong reactions with hFSH and beta hFSH, but no cross reactivity with hTSH, hLH, hCG, alpha hCG, alpha hFSH. The recognition site is located near the alpha 1 binding site on the hFSH molecule. It is not compatible with other anti hFSH beta antibodies.
经测试应用	适用于: IHC-Fr, ELISA, IHC-P
种属反应性	与反应: Human
免疫原	Full length native protein (purified) (Human).
表位	Directed against the Beta 2 epitope of the Beta subunit of human FSH beta.
阳性对照	Pituitary Gland.
常规说明	<p>This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆

克隆编号	INN-hFSH-60
骨髓瘤	P3x63-Ag8.653
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab11069于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

靶标

相关性 Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis. Members of the pituitary glycoprotein hormone family, of which FSH is one consist of a shared alpha chain and a beta chain encoded by a separate gene. The FSH beta gene encodes the beta subunit of follicle stimulating hormone.

细胞定位 Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors